BioCer


HaemoCer™

Endoscopic Applicator


BioCer

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Fases de Coagulación

BioCer

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
BioCer

Apply
HaemoCer PLUS


55% blood plasma
(of which approx. 90% water and 10% proteins)

2% white blood cells (leukocytes) and platelets (thrombocytes)

43% red blood cells (erythrocytes)

Updates in Surgery (2021) 73:1267–1273
https://doi.org/10.1007/s13304-021-01136-x

POSITION PAPER

An international multidisciplinary peer-driven consensus on the optimal use of hemostatic powders in surgical practice

C. Eden[1] · O. C. Buonomo[2] · J. Busch[3] · J. Gilabert-Estellés[4] · R. Medrano[5] · M. Nosotti[6] · S. Renaud[7] · I. Walker[8] · N. Williams[9] · Massimo Carlini[9]

Received: 13 July 2021 / Accepted: 18 July 2021 / Published online: 28 July 2021
© Italian Society of Surgery (SIC) 2021

Abstract
Bleeding is a consequence of insufficient hemostasis and excessive bleeding at a surgical site is associated with an increased risk of post-operative infection, transfusion and re-operation, in addition to increased hospital length of stay and costs. Surgeons employ a range of methods to achieve hemostasis, including topical hemostatic agents of differing composition and properties. Hemostatic powders are a sub-group of topical hemostats, which can be used in helping as adjuncts to manage troublesome bleeding in a variety of situations. As this technology is relatively new and potentially not well known by the broad surgical community, no specific guidelines or recommendations for the optimal use of hemostatic powders in surgery currently exist. A steering group throughout Europe of multidisciplinary surgeons, expert in hemostasis and hemostatics, identified from literature and from personal experience, five key topics. When to use hemostatic powder, the evidence for use, benefits of use, safety remarks and considerations in various surgical specialties. Thirty-seven statements were subsequently drawn from these five key topics. An online survey was sent to 128 high-volume surgeons working in breast surgery, gynaecological and obstetric surgery, general and emergency surgery, thoracic surgery and urological surgery in Europe to assess agreement (consensus) with these statements. Consensus was defined as high if ≥ 75% and very high if ≥ 90% of respondents agreed with a statement. A total of 79 responses were received and consensus among the surgical experts was very high in 27 (73%) statements, high in 8 (22%) statements and was not achieved in 2 (5%) statements. Based on the consensus scores, the steering group produced 16 key recommendations which they considered could improve patient outcomes by reducing post-operative bleeding and its associated complications using hemostatic powder.

Keywords Surgery · Bleeding · Hemostasis · Topic hemostatics

Introduction

Bleeding is an inevitable consequence of surgery and in most cases can be controlled without any adverse consequence to serious morbidity [1–3] and in extreme cases to mortality [4].

Preoperative prevention includes withdrawal of drugs increasing the bleeding risk prior to surgery, including non-

[4] Department of Paediatrics, Obstetrics and Gynecology, Hospital General Universitario Valencia, Valencia, Spain

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
Eden et al. (2021); An international multidisciplinary peer-driven consensus on the optimal use of haemostatic powders in surgical practice. Updates in Surgery. https://doi.org/10.1007/s13304-021-01136-x

Location of GI-bleeding

BioCer


Upper GI	80-90%
Middle GI	5%
Lower GI	10-20%


Frequency and localisation of gastrointestinal bleeding (modified according to Block et al. 2005)

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
n = 2,432,088

US nationwide analysis

Wuerth et al. (2017) Changing Epidemiology of Upper Gastrointestinal Hemorrhage in the Last Decade: A Nationwide Analysis. Springer; Digestive Diseases and Sciences https://doi.org/10.1007/s10620-017-4882-6

Bleeding severity of non-variceal ulcers

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ESGE European Society of Gastrointestinal Endoscopy
Working together in endoscopy to achieve optimal patient care

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Guidelines refer to Forrest classification

ESGE Recommendations

BioCer

Stage	Characteristics	
Forrest Ia	Spurting hemorrhage	Intervention necessary
Forrest Ib	Oozing hemorrhage	
Forrest IIa	Visible vessel	
Forrest IIb	Adherent clot	Remove clot, then re-assess
Forrest IIc	Hematin on ulcer base	No intervention necessary
Forrest III	Clean ulcer base	

HaemoCer Plus GI bleeding algorithm

BioCer

NVUGIB

LGIB

POST INT

Forrest classification	1st-line	adjuvant
Ia	X	✓
Ib	✓	✓
IIa	✓	✓

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
Summary bleeding algorythm

BioCer

Bleeding cohort	Primary bleeding	Post-interventional bleeding
Bleeding Severity	Forrest (Ia), Ib, IIa	
Indications	Peptic ulcers MV syndrome Haemorrhagic oesophagitis Tumours	EMR ESD Polypectomy Biopsy

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

Endoscopic Set

3 Components

1. Application Catheter with Filter
 transports HaemoCer PLUS through the instrument channel
2. Pump Ball
 generates the air flow
3. HaemoCer PLUS
 well-established topical hemostat

BioCer

Manual air-flow generation for powder transport

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

BioCer

Compatible with standard endoscopes

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Preparation and application

BioCer

Connect the pump ball to the Luer Lock connection at the filter end of the applicator.

Use the cap to close the attachment on the Y-Piece.

While the applicator is being inserted, the pump ball should also be used occasionally.

Once the applicator is positioned, keep pumping the pump ball 8-10 times quickly and firmly.

Remove the cap and bring the applicator in an upright position, so that the opening of the bellows bottle faces upwards when it is screwed on.

Turn the applicator by 180° and begin administering the powder by pumping the pump ball several times.

Overview SKU

BioCer

Applicator sets

Item number	Description	Contents/PU
HFS323	HaemoCer PLUS Endoscopic Set 3g	1 Standard Applicator HaemoCer PLUS 3g 1 Endoscopic Applicator